## bioresorbable matrix barrier

Contains One GUIDOR Matrix  
Contenu : Une membrane GUIDOR  
Contiene una barrera matriz GUIDOR

### Rx Only

For Professional Use Only  
Usage professionnel uniquement  
Solo para uso profesional

### STERILE

Sterile Product (e-beam) Do Not Resterilize  
Produit stérile (e-beam) Ne pas restériliser  
Producto esterilizado por haz de electrones. No reesterilizar.

For Single Use Only  
Pour usage unique seulement  
Solo para dosis única

36°F / 2°C 46°F / 8°C

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### English

## INDICATIONS FOR USE

To aid in bone regeneration and augmentation in oral surgery for:

- Extraction socket site preservation
- Immediate implant placement at time of extraction or delayed placement when additional bone regeneration is desired
- Ridge augmentation
- Sinus elevation
- Stable barrier for the containment of bone grafting materials

## CONTRAINDICATIONS

The GUIDOR Matrix Barrier is contraindicated in those situations where:

- General oral surgery should not be performed,
- There is any known sensitivity to these materials (see MATERIAL DESCRIPTION) and/or
- There is active or latent infection at or around the surgical site.

Currently, there are no known additional contraindications to the use of the GUIDOR Matrix Barrier.

## PRECAUTIONS

- The GUIDOR Matrix Barrier is not intended for use in defects other than those stated under INDICATIONS FOR USE.
- The GUIDOR Matrix Barrier has not been clinically tested:
  - For extensive bone augmentation
  - For use in the treatment of failing implants
- In patients with any systemic disorder or disease that involves an unacceptable increase in post operative risk for complications
- In pediatric patients, pregnant or nursing women.
- The GUIDOR Matrix Barrier cannot be reused or resterilized. Resterilization will damage the GUIDOR Matrix Barrier.
- Do not use if foil pouch has previously been opened or damaged.

## ADVERSE REACTIONS

Possible complications following any oral surgery include thermal sensitivity, flap sloughing, some loss of crestal bone height, abscess formation, infection, pain, swelling, and complications associated with the use of anesthesia.

As with any type of surgical therapy, the patient may experience discomfort for a few days.

## MATERIAL DESCRIPTION

The GUIDOR Matrix Barrier is made of a translucent bioresorbable material consisting of poly-D, L-lactide, and poly-L-lactide, along with acetyl-tributylcitrate (ATBC), a citric acid ester.

At storage temperatures 36°F to 46°F (2°C to 8°C), the GUIDOR Matrix Barrier is stiff and brittle and may fracture if manipulated prior to reaching room temperature.

## DESIGN DESCRIPTION


Figure 1. Profile of the GUIDOR Matrix Barrier.

The GUIDOR Matrix Barrier, as shown in Figure 1, consists of a multi-layer construct:

- An external layer (A) which is placed to face the gingival tissue.
- An internal layer (B) which is placed to face the bone.
- The two layers are separated by spacers (C) to form an inner space into which tissue can migrate.

The external layer has rectangular perforations (D) (400 – 500 per cm<sup>2</sup>) that allow gingival tissue to penetrate the inner space for GUIDOR Matrix Barrier stabilization. The internal layer has circular perforations (E) (4000 – 5000 per cm<sup>2</sup>) that act as a barrier to prevent penetration of the gingival tissue until the GUIDOR Matrix Barrier is resorbed and allows tissue integration.

The internal layer has outer spacers (F) that maintain space for bone regeneration.

## PACKAGING

The GUIDOR Matrix Barrier is provided sterile (e-beam) in a transparent plastic tray designed to protect the GUIDOR Matrix Barrier from physical damage. The plastic tray is packed in a sealed foil pouch, which maintains the sterility of the GUIDOR Matrix Barrier. Do not use if the foil pouch has previously been opened or damaged. A variety of sizes are available to meet various clinical needs in the treatment of bone defects.

## WARNINGS

The GUIDOR Matrix Barrier is sterilized (e-beam) and intended for single use only. Do not resterilize or reuse. The GUIDOR Matrix Barrier is not intended for immediate load-bearing.

Do not compromise blood supply to the defect area. Do not use if foil pouch has previously been opened or damaged, or if the expiration date has been exceeded.

## INSTRUCTIONS FOR USE

### 1. PREPARATION OF GUIDOR MATRIX BARRIER PRIOR TO USE

Allow GUIDOR Matrix Barrier to reach room temperature prior to use. Remove the GUIDOR Matrix Barrier from the refrigerator.

- Remove foil pouch from the carton. The Matrix Barrier will reach room temperature in the foil pouch after 15 minutes.
- OR**
- Remove the transparent plastic tray containing the Matrix Barrier from the foil pouch. The GUIDOR Matrix Barrier will reach room temperature in the plastic tray (outside the foil pouch) after 1 minute.

Do not attempt to warm the GUIDOR Matrix Barrier in your hands.

At body temperature, the GUIDOR Matrix Barrier becomes malleable within seconds. The GUIDOR Matrix Barrier should only be placed onto a sterile field. Aseptic handling techniques are essential during application.

## 2. SPACE MAKING

The GUIDOR Matrix Barrier is malleable at body temperature which facilitates application and adaptation of the Matrix Barrier to the adjoining bone surface. Without support i.e., sufficient amount of bone grafting material, the GUIDOR Matrix Barrier may collapse into the defect. The use of appropriate space making support for the Matrix Barrier is left to the judgment of the dental practitioner.

### 3. PREPARATION OF SITE

- It is desirable to use a flap management technique that ensures total defect coverage, if possible.
- Remove all granulation tissue from the bone surface.
- When placing bone grafting materials, follow the manufacturer's instructions for preparation and insertion.

### 4. PLACEMENT OF THE GUIDOR MATRIX BARRIER


Figure 2. Top view of the GUIDOR Matrix Barrier

- Maintain aseptic handling throughout the procedure.
- GUIDOR Matrix Barrier is packaged in the plastic tray with the external layer facing the opening lid. As shown in Figure 2, the external layer is marked with an “S” and multiple 2 mm square notches along the edge of the Matrix Barrier.
- Remove the GUIDOR Matrix Barrier from the plastic tray using a flat tweezer.
- To shape, hold the GUIDOR Matrix Barrier gently between the fingers and cut with sharp scissors. Make sure the Matrix Barrier has reached room temperature before shaping.
- When shaping the GUIDOR Matrix Barrier, the “S” marker and notches may get cut off. Take note of the direction and location of the “S” marker and notches prior to shaping.
- Place grafting material for appropriate space making support.

- Place the GUIDOR Matrix Barrier over the defect area. The external layer with the “S” mark should face the gingiva. Ensure that a sufficient amount of the Matrix Barrier extends over the bone surrounding the defect.
- Avoid applying excessive pressure on the GUIDOR Matrix Barrier.

### 5. FLAP MANAGEMENT

- Securely tuck the Matrix Barrier under the flaps. Whenever possible, the soft tissue flaps should be completely closed and sutured securely but without tension.
- There may be cases where complete closure may not be achievable, e.g. for wide socket openings. In cases where primary closure is not achievable, the GUIDOR Matrix Barrier may be left exposed, provided the Matrix Barrier is securely tucked under the flap and sutures are placed over the exposed Matrix Barrier to secure it and facilitate positioning.<sup>1</sup>

### 6. POST-SURGICAL CARE

#### Instructions for all surgically treated area(s):

- Following surgery, the patient should refrain from mechanical cleaning techniques at or close to the treated area for a minimum of 4 weeks after surgery or until the GUIDOR Matrix Barrier no longer remains visible. During this period, mouth rinsing with an antimicrobial agent, preferably alcohol free, such as chlorhexidine, is recommended.
- Avoid pressure by palpitation or a prosthesis on the treated area during the healing period.
- Antibiotic therapy may be provided at the discretion of the practitioner. If infection or abscess develops, removal of the GUIDOR Matrix Barrier may be indicated depending on the severity of the complication.

#### In cases where the GUIDOR Matrix becomes unintentionally or is intentionally exposed:

- More frequent patient follow-up may be indicated.
- Topical use of an antimicrobial agent, preferably alcohol free, such as chlorhexidine, is strongly recommended for at least 4 weeks or until the GUIDOR Matrix Barrier is no longer present. Topical use of an antimicrobial agent may cause discoloration (darkening) of the Matrix Barrier.
- The use of a systemic antibiotic may be recommended.

- During healing, the GUIDOR Matrix Barrier may fold upon itself. This is expected and does not require removal or trimming of the Matrix Barrier. Trimming of the Matrix Barrier is not advised, unless it relieves patient discomfort. Trimming the membrane may promote early dislodgement causing diminished regenerative results.
- To avoid premature degradation of the GUIDOR Matrix Barrier, the patient should refrain from mechanical cleaning techniques and chewing at or close to the treated area.
- Following exposure, the GUIDOR Matrix Barrier tends to resorb more quickly (in approximately 4 weeks) than if fully submerged.

### RE-ENTRY

The area(s) treated for bone regeneration or augmentation should not be uncovered for approximately 4 to 6 months after placement of the GUIDOR Matrix Barrier. Remnants of the Matrix Barrier may be present at the time of re-entry, usually within the soft tissue flap. These remnants typically have a liquefied and/or granular appearance. No special measures need to be taken, since the remnants will resorb and disappear.

<sup>1</sup>. Rosen PS, Rosen AD. Purposeful Exposure of a Poly(lactic Acid) Barrier to Achieve Socket Preservation for Placement of Dental Implants: Case Series Report. *Compendium* 2013; Vol. 34, No.1: 2-7.

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The GUIDOR name and logo are registered trademarks of Sunstar Suisse SA / Le nom et le logo GUIDOR sont des marques déposées de Sunstar Suisse SA / El nombre y el logotipo de GUIDOR son marcas registradas de Sunstar Suisse SA