SUNSTAR

GUIDOR

bioresorbable matrix barrier

Contains One GUIDOR Matrix
Contenu : Une membrane GUIDOR
Contiene una barrera matriz GUIDOR

Rx Only

For Professional Use Only
Usage professionnel uniquement
Solo para uso profesional

STERILE B

Sterile Product (e-beam) Do Not Resterilize
Produit stérile (e-beam) Ne pas restériliser
Producto esterilizado por haz de electrones. No reesterilizar.

For Single Use Only
Pour usage unique seulement
Solo para dosis única

36°F / 2°C

46°F / 8°C

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English

INDICATIONS FOR USE

The GUIDOR® Bioresorbable Matrix Barrier is indicated for the surgical treatment of periodontal defects, to aid in the regeneration and integration of periodontal tissue components. For use as an adjunct in periodontal surgical treatment to supplement the reparative process following scaling and root planing for:

- Class II furcations
- Intraony defects
- Recession type defects

CONTRAINDICATIONS

The GUIDOR Matrix Barrier is contraindicated in those situations where:

- General oral surgery should not be performed,
- There is any known sensitivity to these materials (see MATERIAL DESCRIPTION), and/or
- There is active or latent infection at or around the surgical site.

Currently, there are no known additional contraindications to the use of the GUIDOR Matrix Barrier.

PRECAUTIONS

- The GUIDOR Matrix Barrier is not intended for use in defects other than those stated under INDICATIONS FOR USE.

- The GUIDOR Matrix Barrier has not been clinically tested:

- For use in the treatment of failing implants
- In patients with any systemic disorder or disease that involves an unacceptable increase in postoperative risk for complications.
- In pediatric patients, pregnant or nursing women.
- In patients with conditions involving extremely severe defects with too little remaining periodontium.

- The GUIDOR Matrix Barrier cannot be reused or resterilized. Resterilization will damage the GUIDOR Matrix Barrier.

- Do not use if foil pouch has previously been opened or damaged.

ADVERSE REACTIONS

Possible complications following any oral surgery include thermal sensitivity, gingival recession, flap sloughing, resorption or ankylosis of the treated root, some loss of crestal bone height, abscess formation, infection, pain, swelling, gingival irregularities, and complications associated with the use of anesthesia. As with any type of surgical therapy, the patient may experience minor discomfort for a few days.

MATERIAL DESCRIPTION

The GUIDOR Matrix Barrier is made of a translucent bioresorbable material consisting of poly-D, L-lactide, and poly-L-lactide, along with acetyl-tributylcitrate (ATBC), a citric acid ester. At storage temperatures 36°F to 46°F (2°C to 8°C), the GUIDOR Matrix Barrier is stiff and brittle and may fracture if manipulated prior to reaching room temperature. It reaches room temperature after 15 minutes in the foil pouch.

At body temperature, the material becomes malleable within seconds.

DESIGN DESCRIPTION


Figure 1. Profile of the GUIDOR Matrix Barrier.

The GUIDOR Matrix Barrier, as shown in Figure 1, consists of a multi-layer construct:

- An external layer (A) which is placed to face the gingival tissue.
- An internal layer (B) which is placed to face the tooth.
- The two layers are separated by spacers (C) to form an inner space into which tissue can migrate. The external layer has rectangular perforations (D) (400 – 500 per cm²) that allow gingival tissue to penetrate the inner space for GUIDOR Matrix Barrier stabilization.

The internal layer has circular perforations (E) (4000 – 5000 per cm²). The internal layer includes outer spacers (F) on the root side of the GUIDOR Matrix Barrier. In transgingival configurations, the GUIDOR Matrix

Barrier has a coronal portion where a bioresorbable ligature (G) is placed. The ligature is used to anchor the GUIDOR Matrix Barrier to the tooth. The top of the internal layer includes a bar (H) which functions as a seal between the GUIDOR Matrix Barrier and the tooth.

PACKAGING

The GUIDOR Matrix Barrier is provided sterile (e-beam) in a transparent plastic tray designed to protect the GUIDOR Matrix Barrier from physical damage. The plastic tray is packed in a sealed foil pouch, which maintains the sterility of the GUIDOR Matrix Barrier. Do not use if foil pouch has previously been opened or damaged.

A variety of straight and curved configurations and sizes are available to meet various clinical needs in the treatment of periodontal defects.

WARNINGS

The GUIDOR Matrix Barrier is sterilized (e-beam) and intended for single use only. Do not resterilize or reuse. Do not compromise blood supply to the defect area. Do not use if foil pouch has previously been opened or damaged, or if the expiration date has been exceeded.

INSTRUCTIONS FOR USE

1. PREPARATION OF GUIDOR MATRIX BARRIER PRIOR TO USE

Allow GUIDOR Matrix Barrier to reach room temperature prior to use. Remove the GUIDOR Matrix Barrier from the refrigerator.

- Remove foil pouch from the carton. The GUIDOR Matrix Barrier will reach room temperature in the foil pouch after 15 minutes.
- OR**
- Remove the transparent plastic tray containing the Matrix Barrier from the foil pouch. The GUIDOR Matrix Barrier will reach room temperature in the plastic tray (outside the foil pouch) after 1 minute.

Do not attempt to warm the GUIDOR Matrix Barrier in your hands. At body temperature, the GUIDOR Matrix Barrier becomes malleable within seconds.

The GUIDOR Matrix Barrier should only be placed onto a sterile field. Aseptic handling techniques are essential during application.

2. PREPARATION OF SITE

- Perform intracrevicular incisions to allow for full thickness flaps.
- Preserve the interproximal papillae in order to ensure future interproximal coverage of the GUIDOR Matrix Barrier.
- Excise pocket epithelium to allow for connective tissue integration with the external layer of the GUIDOR Matrix Barrier.
- Remove all granulation tissue.
- Scale and plane the root surfaces.
- Use of grafting material is at the discretion of the practitioner based on the defect morphology.

3. PLACEMENT OF THE GUIDOR MATRIX

- GUIDOR Matrix Barrier is packaged in the plastic tray with the external layer facing the opening lid. As shown in Figure 2, the external layer is marked with an “S” and multiple 2 mm square notches along the edge of the Matrix Barrier.


Figure 2. Top view of the GUIDOR Matrix Barrier

- Maintain aseptic handling throughout the procedure.
- Make sure that the GUIDOR Matrix Barrier has reached room temperature before manipulation or trimming. Do not attempt to warm the GUIDOR Matrix Barrier in your hands.

- Identify the external layer on the GUIDOR Matrix Barrier. In the plastic tray the external layer faces the lid. The external layer can be identified by the “S” marked on this side, as shown in Figure 2, and may also be marked by a 2 mm square notch(s) in the margin of this layer.

- After the GUIDOR Matrix Barrier reaches room temperature, remove it from the foil pouch and open plastic tray using a flat tweezer.

- If folding of the GUIDOR Matrix Barrier for interproximal insertion is needed, ensure that the GUIDOR Matrix Barrier has obtained room temperature prior to folding.

- To shape, hold the GUIDOR Matrix Barrier gently between the fingers and cut with sharp scissors. Make sure the GUIDOR Matrix Barrier has reached room temperature before shaping. When shaping the GUIDOR Matrix Barrier, the “S” marker and notches may get cut off. Take note of the direction and location of the “S” marker and notches prior to shaping.
- Place the GUIDOR Matrix Barrier to a depth sufficient to cover the defect and a sufficient area of the surrounding tissue to avoid GUIDOR Matrix Barrier collapse into the defect.

- The external layer with the “S” mark should face the gingiva and be placed with the “S” mark facing in the apical direction.
- On application, tie the ligature around the root of the tooth and secure it with a square knot. Do not make an extra sling as this may unlock the knot.

- Overextension of the GUIDOR Matrix Barrier is not expected to have adverse effects.
- Cut the ligature endings. Longer endings can be easier to place under the flap.

- The flaps are sutured to complete coverage of the GUIDOR Matrix Barrier. The choice of suture, slowly resorbing or non-resorbing, is left to the discretion of the practitioner.

4. WOUND CLOSURE

- Make every attempt to achieve primary closure to promote primary healing.
- Avoid tension in the repositioned flap.

5. GUIDOR MATRIX EXPOSURE

- Exposed GUIDOR Matrix Barrier material does not need to be trimmed. It will disappear 6 to 8 weeks after surgery due to combined resorption and attrition.
- The patient should be scheduled for post-surgical follow-up at least every two weeks during the initial 6 weeks of healing.
- After the initial 6 weeks of healing, follow-up visits should be scheduled to supervise and monitor the patient’s plaque control.

- Following exposure, the GUIDOR Matrix Barrier tends to resorb more quickly than if fully submerged.
- To minimize the risk of bacterial infection, application of an antimicrobial agent, preferably alcohol free, such as chlorhexidine, is recommended until the GUIDOR Matrix Barrier is no longer present. Topical use of an antimicrobial agent may cause discoloration (darkening) of the Matrix Barrier.
- To avoid premature degradation of the GUIDOR Matrix Barrier, the patient should refrain from mechanical cleaning techniques and chewing at or close to the treated area.

6. USE OF ANTIBIOTICS

Antibiotic therapy may be provided at the discretion of the practitioner and should follow current standards of care. In cases of infection or abscess, removal of the GUIDOR Matrix Barrier may be indicated depending on the severity of the complication.

7. POST-SURGICAL CARE

Special instructions are needed for the surgically treated area(s). It is of utmost importance for the success of the procedure, that the area is kept free from contamination. Care should also be taken not to disrupt the wound by mechanical means.

- Following surgery the patient should refrain from mechanical cleaning techniques near the treated area for a minimum of 4 weeks after surgery or until the GUIDOR Matrix Barrier no longer remains visible. During this period, mouth rinsing with an antimicrobial agent, preferably alcohol free, such as chlorhexidine, is recommended.

- The patient should be scheduled for post-surgical follow-up at least every two weeks during the initial 6 weeks of healing.
- After the initial 6 weeks of healing, follow-up visits should be scheduled to supervise and monitor the patient’s plaque control.

- Supragingival professional prophylaxis can be performed as needed. However, do not use pumice or other types of polishing paste during the initial 4 weeks of healing.
- Subgingival scaling should not be performed until at least 6 months after surgery.
- Probing the defect for treatment evaluation should not be performed until at least 6 months after surgery. Retreatment may be considered.

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The GUIDOR name and logo are registered trademarks of Sunstar Suisse SA / Le nom et le logo GUIDOR sont des marques déposées de Sunstar Suisse SA / El nombre y el logotipo de GUIDOR son marcas registradas de Sunstar Suisse SA